


# The QLK-100

## Bracket Reference List

This bracket reference list is meant to be used as a general application guide, and does not guarantee that a bracket will be available for your wheelchair. The bracket design and development is based on the specifications obtained from manufacturers on the basic model of each wheelchair, and it does not take into consideration any additional accessories or customized features.



Please have the following information available when visiting your QLK-100 dealer:

- Manufacturer of the wheelchair
- Model or serial number
- Take note of any special features or accessories (i.e. tilt, seat modifications, elevating power seats, etc.) as these may affect bracket installation.
- Confirm tubing diameter in manual wheelchair applications.
  - Minimum ground clearance required is 3 inches beneath the wheelchair.
  - Driver application require front stabilizer base.
  - Wheelchairs with backrests below the shoulder blades are not recommended to be used with the QLK-100.
- Confirm length of QLK wire harness required for application.

**Wire harness available in the following lengths:**

**Q04C0138** -20" longer (115" in total length) \$18.00

**Q04C0083**-95" long, standard length included with the QLK 100 \$18.00

## Manual Wheelchairs

### Invacare

Model Name	Specifications	Q-Bracket	Reference
9000 Series	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
9000 XT	Single frame	Q04BR001	BKT-Q2
9000 XT	Bi-level frame	Q04BR003	BKT-XT
Action Allegro	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Action A-T	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Action Comet	Frame rail less 6" from floor	Q04BR003	BKT-9000-1
Action Comet	Frame rail above 6" from floor	Q04BR001	BKT-Q2
Action Junior	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Action MVP	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Action Orbit	Frame rail less 6" from floor	Q04BR003	BKT-9000-1
Action Orbit	Frame rail above 6" from floor	Q04BR001	BKT-Q2
Action Patriot	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Action Super Pro-T	Frame rail less 6" from floor	Q04BR003	BKT-9000-1
Action Super Pro-T	Frame rail above 6" from floor	Q04BR001	BKT-Q2
Action Tiger	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Action Xtra	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
ProSPIN X4		Q04BR001	
Ride Lite 9000	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Rolls 2000	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Rolls 4000	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Rolls 900	Frame rail less 6" from floor	Q04BR004	BKT-9000
Rolls 900	Frame rail above 6" from floor	Q04BR001	BKT-Q278
Solara 2nd generation	May have to modify bracket depending on w/c specs; If chair has less than 5" in ground clearance	Q04BR006	BKT-SOLO-2G
Solara 2nd generation	May have to modify bracket depending on w/c specs; If chair has more than 5" in ground clearance. Customer will also need 4 pieces of Q04C0036	Q04BR001	
Solo / Solara Jr./Spree GT	1" tubing diameter/Frame rail less 6" from floor	Q04BR003	BKT-9000-1
Spree		Q04BR003	BKT-9000-1
Tracer Series	Frame rail less 6" from floor (7/8" Tubing)	Q04BR002	BKT-9000
Tracer Series	Frame rail above 6" from floor (1" Tubing)	Q04BR003	BKT-Q278
Tracer Series	Frame rail above 6" from floor (7/8" Tubing)	Q04BR004	
Tracer LT		Q04BR002 or 004	N/A
Tracer LX		Q04BR002 or 004	N/A
Youthmobile 3000	1" tubing diameter/ frame rail less 6" from floor	Q04BR003	BKT-9000-1
Youthmobile 9000	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2

### Sunrise Medical (Quickie)

Model Name	Specifications	Q-Bracket	Reference
Breezy	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Breezy 2		Q04BR001	BKT-Q2
GP / GPV	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Iris		Q04BR052	BKT-IRIS
LX Series	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Quickie 2	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2

# The **QLK-100**

## Bracket Reference List

This bracket reference list is meant to be used as a general application guide, and does not guarantee that a bracket will be available for your wheelchair. The bracket design and development is based on the specifications obtained from manufacturers on the basic model of each wheelchair, and it does not take into consideration any additional accessories or customized features.

Quickie 2HP	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Quickie TS	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Quickie TS	1" tubing diameter/ frame rail is at/or below 4" from floor	Q04BR003	BKT-Q2
Recliner	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Zippie	Frame rail less 6" from floor	Q04BR003	BKT-9000-1
Zippie	Frame rail above 6" from floor	Q04BR001	BKT-Q2
Zippie 2		Q04BR001	BKT-Q2
Zippie GS		Q04BR001 or 003	N/A
Zippie TS		Q04BR001	BKT-Q2

# The QLK-1000

## Bracket Reference List

This bracket reference list is meant to be used as a general application guide, and does not guarantee that a bracket will be available for your wheelchair. The bracket design and development is based on the specifications obtained from manufacturers on the basic model of each wheelchair, and it does not take into consideration any additional accessories or customized features.

### Other

Model Name	Specifications	Q-Bracket	Reference
Bantex	7/8" tubing diameter/ frame rail above 6" from floor	Q04BR003	BKT-Q278
Bantex	7/8" tubing diameter/ frame rail less 6" from floor	Q04BR002	BKT-9000
Comfort Orthopedic	7/8" tubing diameter/ frame rail above 6" from floor	Q04BR003	BKT-Q278
Comfort Orthopedic	7/8" tubing diameter/ frame rail less 6" from floor	Q04BR002	BKT-9000
Drive Medical - Cirrus III		Q04BR001	BKT-Q2
Drive Medical - Cirrus IV		Q04BR001	BKT-Q2
Freedom Design - Libre	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Gendron	7/8" tubing diameter/ frame rail above 6" from floor	Q04BR003	BKT-Q278
Gendron	7/8" tubing diameter/ frame rail less 6" from floor	Q04BR002	BKT-9000
Gunnell	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
Jaken Medical	7/8" tubing diameter/ frame rail above 6" from floor	Q04BR002 or 004	BKT-Q278
Jaken Medical	7/8" tubing diameter/ frame rail less 6" from floor	Q04BR002	BKT-9000
Kuschall Champ 1000	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
LaBac Manual	7/8" tubing diameter/ frame rail less 6" from floor	Q04BR002	BKT-9000
Lumex	7/8" tubing diameter/ frame rail above 6" from floor	Q04BR002 or 004	BKT-Q278
Lumex	7/8" tubing diameter/ frame rail less than 6" from floor	Q04BR002	BKT-9000
Merit		Q04BR002	N/A
Merit P310 Series		Q04BR032	
Summit	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
TiLite - Folding	1" tubing diameter/ frame rail above 6" from floor	Q04BR001	BKT-Q2
TiLite Rigid Evo		Q04BR001	
TiLite Folding 2GX	1" tubing diameter/frame rail above 6" from floor	Q04BR001	BKT-Q2
TiLite Folding Aero X	1" tubing diameter/frame rail above 6" from floor	Q04BR001	BKT-Q2
TiLite Folding TX	1" tubing diameter/frame rail above 6" from floor	Q04BR001	BKT-Q2
TiLite Folding SX	1" tubing diameter/frame rail above 6" from floor	Q04BR001	BKT-Q2
TiLite Folding X	1" tubing diameter/frame rail above 6" from floor	Q04BR001	BKT-Q2
TiLite Youth YGS	1" tubing diameter/frame rail above 6" from floor	Q04BR001	
TiLite Youth YGS	1" tubing diameter/frame rail less than 6" from floor	Q04BR003	
Medline Excel 1000	1" tubing diameter/frame rail less than 6" from floor	Q04BR003	
Medline Excel 2000	1" tubing diameter/frame rail less than 6" from floor	Q04BR003	
Medline Excel K1	1" tubing diameter/frame rail less than 6" from floor	Q04BR003	
Medline Excel K4	1" tubing diameter/frame rail above 6" from floor	Q04BR001	

www.qstraint.com

# The **OLK-1000**

## Power Wheelchairs

### Invacare

Model Name	Specifications	Q-Bracket	Reference
Arrow (2nd Generation)	Std. motors & model starts with 2G, rear & mid drive	Q04BR008	BKT-2G
Arrow (2nd Generation)	Gearless brushless motors & model starts with 2G	Q04BR009	BKT-GB2
Arrow (3rd Generation)	Std. motors & model starts with 3G, rear & mid drive	Q04BR008	BKT-2G
Arrow (3rd Generation)	Gearless brushless motors & model starts with 3G	Q04BR010	BKT-GB3
M51	Forward entry only	Q04BR018	BKT-M61
M61	Forward entry only	Q04BR018	BKT-M61
M71	Forward entry only	Q04BR018	BKT-M71
M91	Forward entry only	Q04BR019	BKT-M91
M94	Forward entry only	Q04BR019	BKT-M91
Pronto M51/M61 (Forward entry only)	Specify foot rest options	Q04BR018	BKT-M61
Pronto M71 (forward entry only)	Specify foot rest options	Q04BR018	BKT-M71
Pronto M91/M94 (Forward entry only)	Specify foot rest options	Q04BR019	BKT-M91
Ranger X (2nd Generation Storm Series)	Stdard 4 pole motors & model starts with 2G	Q04BR008	BKT-2G
Ranger X (2nd Generation Storm Series)	Gearless brushless motors & model starts with 2G	Q04BR009	BKT-GB2
Ranger X (3rd Generation Storm Series)	Std. 4 pole motors & model starts with 3G	Q04BR008	BKT-2G
Ranger X (3rd Generation Storm Series)	Gearless brushless motors & model starts with 3G	Q04BR010	BKT-GB3
TDX 3 (Forward entry only)	2 group 22 batteries&standard seating	Q04BR071	BKT-TDX222
TDX 4 (Forward entry only)	2 group 22 batteries&standard seating	Q04BR071	BKT-GB4
TDX 5 (Forward entry only)	2 group 22 batteries&standard seating	Q04BR071	BKT-TDX223
TDX 3 (Forward entry only)	2 group 24 batteries&standard seating	Q04BR070	BKT-TDX224
TDX 4 (Forward entry only)	2 group 24 batteries&standard seating	Q04BR070	BKT-TDX224
TDX 5 (Forward entry only)	2 group 24 batteries&standard seating	Q04BR070	BKT-TDX224
TDX 3 (Forward entry only)	Motion Concept Seat, 2 group 22 batteries	Q04BR060	BKT-TDX-MC
TDX 4 (Forward entry only)	Motion Concept Seat, 2 group 22 batteries	Q04BR060	BKT-2G
TDX 5 (Forward entry only)	Motion Concept Seat, 2 group 22 batteries	Q04BR060	BKT-TDX-MC
TDX 3 (Forward entry only)	Invacare Formula Invisible Super Low Tilt seating	Q04BR068	BKT-TDX-SLT
TDX 4 (Forward entry only)	Invacare Formula Invisible Super Low Tilt seating	Q04BR068	BKT-TDX-SLT
TDX 5 (Forward entry only)	Invacare Formula Invisible Super Low Tilt seating	Q04BR068	BKT-TDX-SLT
TDX 3 (Forward entry only)	TRE Seating System	Q04BR069	BKT-TDX-TRE
TDX 4 (Forward entry only)	TRE Seating System	Q04BR069	BKT-TDX-TRE
TDX 5 (Forward entry only)	TRE Seating System	Q04BR069	BKT-TDX-TRE
TDX SP	Group 22 Batteries (Works with or without power seat) Group 22 battery base has a width of 24 inches.	Q04BR058	
TDX SP	Group 24 Batteries (Group 24 battery base has a width of 25.5 inches )	Q04BR059	
Torque (2nd Generation w/standard 4 pole motors)	Rear & mid wheel drive. Model starts 2G	Q04BR008	BKT-2G
Torque (2nd Generation w/gearless brushless mtr)	Rear & mid wheel drive. Model starts 2G	Q04BR009	BKT-GB2
Torque (3rd Generation w/standard 4 pole motors)	Rear & mid wheel drive. Model starts 3G	Q04BR008	BKT-2G
Torque (3rd Generation w/gearless brushless mtr)	Rear & mid wheel drive. Model starts 3G	Q04BR010	BKT-GB3

### Sunrise Medical Quickie

Model Name	Specifications	Q-Bracket	Reference
Freestyle F11 (Forward entry)	Not to be equipped with the high back comfort seat	Q04BR037	BKT-F11
Freestyle M11 (Forward entry)	Not to be equipped with the high back comfort seat	Q04BR038	BKT-M11
P200	Power Configuration	Q04BR022	BKT-QP2
P210	Power Configuration	Q04BR022	BKT-QP2
P220	Power Configuration	Q04BR044	BKT-QP222
P222	Power Configuration	Q04BR044	BKT-QP222
P300	20" rear wheels. Round frame tubing	Q04BR003	BKT-9000-1
Rhythm	Forward entry suggested	Q04BR057	BKT-RTM
Groove	Front wheel drive and Rear wheel drive	Q04BR056	
S-525		Q04BR049	BKT-525
S-626	Not to be equipped with elevating seat	Q04BR050	BKT-626
S-646	Not to be equipped with elevating seat	Q04BR050	BKT-626

### Pride Mobility

Model Name	Specifications	Q-Bracket	Reference
Blast 650	May use power elevating seat	Q04BR035	BKT-BLAST
Blast 850	May use power elevating seat	Q04BR035	BKT-BLAST
Blast HD	May use power elevating seat	Q04BR035	BKT-BLAST
Jazzy 600	May use power elevating seat	Q04BR028	BKT-J600
Jazzy 1100	May use power elevating seat	Q04BR032	BKT-JAZZY
Jazzy 1101	May use power elevating seat	Q04BR033	BKT-J2101
Jazzy 1120	May use power elevating seat	Q04BR032	BKT-JAZZY
Jazzy 1121	May use power elevating seat	Q04BR033	BKT-J2101
Jazzy 1122	May use power elevating seat	Q04BR032	BKT-JAZZY
Jazzy 1143	May use power elevating seat	Q04BR034	BKT-J1143
Jazzy 1143 Ultra	May use power elevating seat	Q04BR028	BKT-J1143-U
Jazzy 1170	May use power elevating seat	Q04BR032	BKT-JAZZY

# The QLK-1000

## Power Wheelchairs

Jazzy 1400	May use power elevating seat	Q04BR032	BKT-JAZZY
Jazzy 1420	May use power elevating seat	Q04BR032	BKT-JAZZY
Jazzy 1470	May use power elevating seat	Q04BR032	BKT-JAZZY
Jazzy 6000	May use power elevating seat	Q04BR028	BKT-J6K
Quantum 600	May use power elevating seat	Q04BR028	BKT-J600
Quantum 6000 & Q6000Z	May NOT use power elevating seat	Q04BR028	BKT-J6K
Quantum 6000 & Q6000Z	May use power Tilt and Recline	Q04BR029	
Quantum 6000T	May use power elevating seat	coming soon	
Quantum Blast 650	May use power elevating seat	Q04BR035	BKT-BLAST
Quantum Blast 850	May use power elevating seat	Q04BR035	BKT-BLAST
Quantum Blast HD	May use power elevating seat	Q04BR035	BKT-BLAST
Quantum Vibe	May use power elevating seat	Q04BR036	BKT-VIBE
Valor	May use power elevating seat	Q04BR036	BKT-VIBE
Vanguard	May use power elevating seat	Q04BR033	BKT-J2101
Vibe R4000	May use power elevating seat	coming soon	

### Permobil

Model Name	Specifications	Q-Bracket	Reference
Chairman Basic	Forward entry	Q04BR029	BKT-CH3
Chairman Basic	Back in Entry	Q04BR029	BKT-CH3-B
Chairman Corpus	Forward entry	Q04BR029	BKT-CH4
Chairman Corpus	Back in Entry	Q04BR029	BKT-CH3-B
Chairman Entra	Forward entry. Must have Permobil clearance package	Q04BR029	BKT-CH3
Chairman Entra	Back in Entry. Must have Permobil clearance package	Q04BR029	BKT-CH3-B
Chairman 2K (Permobil C2K)		Q04BR030	BKT-C2K
Chairman 2S (Permobil C2S)		Q04BR030	BKT-C2K
C300		Q04BR061	BKT-C300
C400		Q04BR030	BKT-C2K
C500		Q04BR032	BKT-C2K
Street		Q04BR013	BKT-STREET
C350		coming soon	

### Other

Model Name	Specifications	Q-Bracket	Reference
Bruno PWC-2200	Fits standard adjustable seat post, PWC-23103	Q04BR012	BKT-PWC-2200
Bruno PWC-2300	Fits standard adjustable seat post, PWC-23103	Q04BR012	BKT-PWC-2300
Teknique FWD - Hoveround	Front wheel drive	Q04BR011	BKT-TEK-FWD
Teknique RWD - Hoveround	Rear wheel drive	Q04BR011	BKT-TEK
Orthofab-Oasis		Q04BR033	
Eagle-Liberty Easy		coming soon	
Electric Mobility P301		coming soon	
Electric Mobility P710		coming soon	
Falcon Rehab Traxx		coming soon	
Golden Technologies Alante		coming soon	
Golden Technologies Compass		coming soon	
IMC Heartway Allure		coming soon	
Innovation in Motion		coming soon	
Levo Combi		coming soon	
Pacesaver Boss 4.5		coming soon	
Pacesaver RFP4		coming soon	
Shoprider Jetstream L		coming soon	
Shoprider Jetstream M		coming soon	
TEFTEC		coming soon	
The Standing Company Power Stander		coming soon	
Tuffcare		coming soon	

www.qstraint.com